FOR ENERGY EFFICIENT INNOVATIONS

**THINK ON.**

[www.onsemi.cn](http://www.onsemi.cn)

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
Figure 2. RSL10 Application Diagram in Buck Mode



Figure 3. RSL10 Application Diagram in LDO Mode

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




Modeler

sound designer modeler map editor library manager control panel calconfig  
TRANSDUCER MODELING

[Create a new .tdr file for the Front Microphone](#)

[Create a new .tdr file for the Rear Microphone](#)

[Create a new .tdr file for the Receiver](#)

[Import a legacy .mlr file and convert to .tdr](#)

[Import a .csv file and convert to .tdr](#)


Double click to open a file from the workspace:

- ▀ Microphones
  - em4346cx.tdr
  - flat\_mic.tdr
- ▀ Receivers
  - ed3146\_10mm.tdr
  - flat\_rec.tdr

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
Create a NEW .tdr file for  
Microphone

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


sound designer modeler map editor library manager control panel calconfig

FITTING PARAMETER MAP EDITOR

Save

Save As

APPLICATION

WIDE DYNAMIC RANGE CO...

FEEDBACK CANCELLER

WIDEBAND AGCO

HARDWARE CONTROL

NOISE REDUCTION

FRONT-END PROCESSING

EQUALIZER

SOUND GENERATOR

MAN MACHINE INTERFACE

DIAGNOSTICS REPORTER

AUTOMATIC RECEIVER DETE...

DATA LOGGER

tone sequencer

FEEDBACK PATH MEASUREMENT...

WIRELESS

FRAMEWORK

ACOUSTIC INDICATORS

ENVIRONMENT CLASSIFIER

Name	Description	Units
X_WDRC_ChannelOutputLimit[n]	The maximum output level for a given channel.	dBFS
X_WDRC_CrossoverFrequency[0]	This is the crossover frequency that separates adjacent channels.	Hz
X_WDRC_CrossoverFrequency[10]	This is the crossover frequency that separates adjacent channels.	Hz
X_WDRC_CrossoverFrequency[11]	This is the crossover frequency that separates adjacent channels.	Hz
X_WDRC_CrossoverFrequency[12]	This is the crossover frequency that separates adjacent channels.	Hz
X_WDRC_CrossoverFrequency[13]	This is the crossover frequency that separates adjacent channels.	Hz
X_WDRC_CrossoverFrequency[14]	This is the crossover frequency that separates adjacent channels.	Hz
X_WDRC_CrossoverFrequency[1]	This is the crossover frequency that separates adjacent channels.	Hz
X_WDRC_CrossoverFrequency[2]	This is the crossover frequency that separates adjacent channels.	Hz

Name: X\_WDRC\_CrossoverFrequency[11] UFN: WDRC Channel Crossover Frequency 11

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



## Control Panel

## Sound Designer 软件之Control Panel 介绍

sound designer modeler map editor library manager control panel calconfig

DEVICE CONTROL PANEL

[Create a new .param file using:](#)

[Create a new .param file by detecting the connected device's Product Library](#)

Double click to open a file from the workspace:

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## Control Panel

sound designer modeler map editor library manager control panel calconfig

DEVICE CONTROL PANEL

Product Library: **Ezairo 7160 SL Demo Library** Product: Ezairo 7160 SL 16 Channels ▾

File

Graphs

Programmer

Memory: Mem A ▾

Noise Reduction	F B C	Framework	Biquads	Wireless	Scratch Memory	Data Log	Hardware	Feedback Measurement	
Application	Diagnostics	S G	AGCo	Front End	Equalizer	E C	MMI	WDRC	Acoustic Indicators
Name		Value		Details					
AUX EQ Enable		Disabled ▾							
Auxiliary Input Crossover Frequency		625		Units: Hz					
Auxiliary Input High Band Gain		0		Units: dB					

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

## Control Panel--- Noise Reduction Option

Application	Diagnostics	S G	AGCo	Front End
Noise Reduction	F B C	Framework	Biquads	Wireless
Name	Value			Details
Noise Reduction Enable	Enabled			 
Noise Reduction Maximum Depth 0	6			Units: dB  
Noise Reduction Maximum Depth 1	6			Units: dB  
Noise Reduction Maximum Depth 2	6			Units: dB  
Noise Reduction Maximum Depth 3	6			Units: dB  
Noise Reduction Maximum Depth 4	6			Units: dB  
Noise Reduction Maximum Depth 5	6			Units: dB  
	6			Units: dB

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## Control Panel--- FBC Option

Application	Diagnostics	S G	AGCo	Front End	E
Noise Reduction	F B C	Framework	Biquads	Wireless	
Name	Value		Details		
FBC Active Sensitivity	13		 		
FBC Active Speed	1		 		
FBC Active Time	1500		Units: milliseconds  		
FBC Enable	Disabled		 		
FBC Gain Management Enable	Disabled		 		
FBC Gain Management Limit	0		Units: dB  		
FBC Idle Speed	7		 		
FBC Startup Time	4		Units: seconds  		

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
## Control Panel--- Wireless Option

Application	Diagnostics	S G	AGCo	Front End
Noise Reduction	F B C	Framework	Biquads	Wireless


Pairing Advertising Mode Timeout	180	Units: seconds
Radio TX Power Level	0 dBm	Units: dBm
Select the binaural connection role	Peripheral	
Slave Latency	0	Units: connection ev
Stream Address	15912422	Range: 0 to 1677721
Stream Key 0	0	Range: 0 to 1677721
Stream Key 1	0	Range: 0 to 1677721
Stream Key 2	0	Range: 0 to 1677721
Stream Key 3	0	Range: 0 to 1677721
Stream Key 4	0	Range: 0 to 1677721

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Calconfig

sound designer modeler map editor library manager control panel calconfig  
CALIBRATION AND CONFIGURATION

[Create a new .clc file using:](#)

Double click to open a file from the workspace:

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## Calconfig

sound designer modeler map editor library manager control panel calconfig

CALIBRATION AND CONFIGURATION

? Edit Programmer Close

Verify the Calibration steps and click Start

Active	Step	Type	Result
<input checked="" type="checkbox"/>	Front Microphone Calibration	Automatic	Pending
<input type="checkbox"/>	Rear Microphone Calibration	Disabled	Skip
<input checked="" type="checkbox"/>	Receiver Calibration	Automatic	Pending
<input type="checkbox"/>	Telecoil Calibration	Disabled	Skip
<input checked="" type="checkbox"/>	Configure Device		Pending
<input checked="" type="checkbox"/>	Burn Scratch Memory		Pending
<input checked="" type="checkbox"/>	Burn Voice Alerts		Pending

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[http://www.jhearing.com/jhear\\_files.htm](http://www.jhearing.com/jhear_files.htm)

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— □ ×

The screenshot shows the JHER software interface. At the top, there is a header bar with the JHER logo on the left and a '连接助听器' (Connect Hearing Aid) button in the center. To the right of the button are two panels for '右耳' (Right Ear) and '左耳' (Left Ear), each displaying '6 memories, 16 Channel ITE/ITC 119' and a 'Memory 0' dropdown menu with a '保存' (Save) button. Below the header is a sidebar menu with icons and labels: '用户管理' (User Management), '助听器型号' (Hearing Aid Model), '自动验配' (Automatic Fitting), '基本验配' (Basic Fitting), '详细验配' (Detailed Fitting), '个性设置' (Personal Settings), '程序设置' (Program Settings), '完成验配' (Complete Fitting), and '验配助手' (Fitting Assistant). The main area is currently blank. At the bottom, there are three status bars: '状态栏1' (Status Bar 1) on the left, '状态栏2:主要显示验配工程中详细的错误信息' (Status Bar 2: Mainly displays detailed error information in the fitting process) in the center, and '当前验配的用户' (Current fitting user) on the right, with the text '您还没有选择用户' (You have not selected a user) below it. A red arrow points from the '连接助听器' button to the text '点这里连接助听器, 读取助听器的参数' (Click here to connect the hearing aid, read the parameters of the hearing aid).

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
# JH-FITTING PC端验配软件 --- RSL10 Dongle 无线验配

# JH-FITTING PC端验配软件 --- RSL10 Dongle 无线验配

The screenshot displays the JH-FITTING PC software interface. At the top, there are controls for connecting a hearing aid, with fields for '右耳' (Right Ear) and '左耳' (Left Ear), each showing '6 memories, 16 Channel ITE/ITC' and a 'Memory 0' dropdown with a '保存' (Save) button.

The main interface is divided into a left sidebar and a main content area. The sidebar contains menu items: 用户管理 (User Management), 助听器型号 (Hearing Aid Model), 自动验配 (Automatic Fitting), 基本验配 (Basic Fitting), 详细验配 (Detailed Fitting), 个性设置 (Personal Settings), 程序设置 (Program Settings), 完成验配 (Complete Fitting), and 验配助手 (Fitting Assistant).

The main content area features a '用户资料' (User Information) section with a search bar, '查找记录' (Search Records), and '新建记录' (New Record) buttons. To the right are buttons for '删除' (Delete), '修改' (Modify), '导入' (Import), and '导出' (Export).

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Overlaid on the table area are four numbered points in red and green text:

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# JH-FITTING PC端验配软件 --- RSL10 Dongle 无线验配

The screenshot displays the JH-FITTING software interface for RSL10 Dongle wireless fitting. It features two frequency response graphs (dB vs. kHz) comparing different fitting levels: 轻声 (Light), 中声 (Medium), and 大声 (Loud). The graphs show the frequency response curves for each level, with the 轻声 level generally having the highest gain and the 大声 level having the lowest gain. The x-axis represents frequency from 125 Hz to 8 kHz, and the y-axis represents gain in dB from 0 to 120.


Below the graphs is the "程序管理" (Program Management) section, which includes a table for managing fitting programs. The table has columns for program ID, name, and actions (Default, Copy, Paste). The table is divided into two sections for "右耳" (Right Ear) and "左耳" (Left Ear).

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Detailed description: This screenshot shows the 'Personal Preferences' (个人喜好) settings screen. It features three sliders for 'Hearing your own voice' (听到自己的声音), 'Sound quality' (声音质量), and 'Sound amplification' (声音放大). Below these is an equalizer graph with three frequency sliders labeled 'Low' (低音), 'Mid' (中音), and 'High' (高音). The graph shows a curve with values of -3dB at the low end, -17dB at the mid end, and -2dB at the high end. A 'Save' (保存) button is located at the top right. The bottom navigation bar contains icons for 'Home' (主页), 'Wireless Link' (无线链接), 'Electronic Medical Record' (电子病历), and 'Remote Assistance' (远程协助).

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Detailed description: This screenshot shows the 'Configuration Information' (配置信息) settings screen. It contains several toggle switches: 'Noise Reduction' (降噪处理) is off, 'Echo Cancellation' (回声处理) is on, 'Wind Noise Switch' (风噪开关) is on, and 'Tinnitus Masking' (耳鸣掩蔽) is on. There is a 'Microphone Settings' (麦克风设置) option with 'Wideband ADM' and a chevron icon. The 'External Audio Switch' (外部音频开关) is also on. At the bottom, there is a 'Volume' (音量) slider. A 'Save' (保存) button is at the top right. The bottom navigation bar is identical to the previous screen, with icons for 'Home' (主页), 'Wireless Link' (无线链接), 'Electronic Medical Record' (电子病历), and 'Remote Assistance' (远程协助).

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